

# BACTERIA TESTS

## DIPSLIDES, TUBE TESTS, INCUBATORS

EASY AND RELIABLE

FLEXIBLE PADDLE



DIPSLIDES FOR VARIOUS SORTS AND COMBINATIONS OF BACTERIA

### VERSIONS

TTC (TOTAL VIALBLE COUNT)

YEAST-MOULDS-FUNGI/TTC

COLIFORMS/TTC

E.COLI - PSEUD. AE./TTC

PSEUDOMON. SPEC./TTC

SULPHATE RED. BACT.

NITRITE RED. BACT/TTC

### BACTEDT057

INCUBATOR FOR 2 DIPSLIDES



### BACTEDT055

INCUBATOR FOR 8 DIPSLIDES



### BACTEDT059

INCUBATOR FOR 10 DIPSLIDES



### BACTEDT052

INCUBATOR FOR 25 DIPSLIDES



### BACTEDT093

INCUBATOR FOR 45 DIPSLIDES



# BACTERIA TESTS

## DIPSLIDES, TUBE TESTS, INCUBATORS

### FEATURES

The definition of a dipslide is a means of testing the microbial content of liquids or surfaces. It consists of a plastic carrier bearing a sterile culture medium which can be dipped in the liquid or pressed on to the surface to be sampled. It is then incubated for 48 hours at a specific temperature to allow microbial growth.

The resulting microbial colonies are estimated by reference to a chart. If the treatment programme is effective, the dipslide count will be consistently low ( $<10^4$ ). If a high result is obtained, the test should be repeated and the treatment programme checked.

Water-i.d.<sup>®</sup> dipslides for bacterial testing are provided as 10 dipslides in a box.

Water-i.d.<sup>®</sup> offers a variety of incubators along with the dipslide range.

#### Water-i.d.<sup>®</sup> Dipslides for Bacteria testing



Example:  
Incubator BactEDT055  
for 8 dipslides

### DIPSLIDES

Each Water-i.d.<sup>®</sup> Dipslide pack contains 10 dipslides along with a comparison chart to obtain the test result.

Item code	Product Description
-----------	---------------------

BactD001	TTC/TTC Dipslide Tests for <b>Total Viable Count</b> in waters and surfaces. Flexible paddle for effective surface contact.
BactD002	Malt/TTC Dipslide Tests for <b>Yeasts, Moulds, Fungi and Total Viable Count</b> in waters and surfaces. Flexible paddle for effective surface contact.
BactD003	RBS/TTC Dipslide Tests for <b>Yeasts, Moulds, Fungi and Total Viable Count</b> for food and manufacturing industry. Flexible paddle for effective surface contact.
BactD004	MAC/TTC Dipslide Tests for <b>Coliforms and Total Viable count</b> in process water. Flexible paddle for effective surface contact.
BactD005	Chromogenic/TTC Dipslide Tests for <b>E.coli, Pseudomonas Aeruginosa and Total Viable Count</b> for enumeration of E.coli. Flexible paddle for effective surface contact.
BactD006	PDM/MAC Dipslide Tests for <b>Pseudomonas species and Coliforms</b> for potable water, pools and spas. Flexible paddle for effective surface contact.
BactD007	PDM/TTC Dipslide Tests for <b>Pseudomonas species and Total Viable Count</b> for Closed water systems, pools and spas. Flexible paddle for effective surface contact.
BactD008	SRB Tube Tests for <b>Sulphate reducing bacteria</b> in Closed water systems. Flexible paddle for effective surface contact.
BactD009	NRB Tube Tests for <b>Nitrite reducing bacteria</b> in Closed water systems. Flexible paddle for effective surface contact.

### INCUBATORS

Item code	Product Description	Voltage
BactEDT057	Incubator for <b>2 Dipslide</b> Tests. Preset to 30°C.	110 / 220 V
BactEDT055	Incubator for <b>8 Dipslide</b> Tests. With temperature gauge and adjustable temperature.	12 / 110 / 220 V
BactEDT059	Incubator for <b>10 Dipslide</b> Tests. With temperature gauge and adjustable temperature. Ideal for NRB/SRB tests. Only offered with 240 V.	12 / 110 / 220 V
BactEDT052	Incubator for <b>25 Dipslide</b> Tests. With temperature gauge and adjustable temperature.	12 / 110 / 220 V
BactEDT093	Incubator for <b>45 Dipslide</b> Tests. With temperature gauge and adjustable temperature.	12 / 110 / 220 V

We gladly Private Label our products for you!  
Write to [sales@water-id.com](mailto:sales@water-id.com) to receive your individual quotation!